

**Meas.Nr.** 69458      **Date** Thursday, 24 October 2013 08:33 AM      **Personal Report No.:**  
**Material** Indigo (MTP + Apollo)      30 min P6cl 80ml 1mm ZrO2      **M130265**  
**Description** pre disp. in ethanol + 1min LABORETTE 17; SVM measurement

**Calculation** Fraunhofer      broad      **TradeOff** 1.000000e+002  
**Scans Fine** 100      **Scans Coarse** 0      **Channels** 51      **Beam Obscuration** 17.00 %  
**Meas. Range** 0.08 [µm] - 42.30 [µm]      **Pump** 0      **Ultrasonics** 0

%	< µm
5	0.41
10	0.52
15	0.59
20	0.66
25	0.72
30	0.78
40	0.90
50	1.02
60	1.16
70	1.33
80	1.53
90	1.87
95	2.16
98	2.48
99	2.72

%	< µm
0.7	0.1
2.3	0.2
9.1	0.5
15.8	0.6
31.9	0.8
48.3	1.0
62.6	1.2
78.5	1.5
88.4	1.8
92.7	2.0
98.2	2.5
99.6	3.0
100.0	3.5
100.0	4.0
100.0	5.0

