

Meas.Nr. 69526 **Date** Monday, 28 October 2013 12:59 PM **Personal Report No.:**
Material Indigo (MTP + Apollo) 45 min P7pl 0,1 mm ZrO2 **M130265**
Description pre disp. in ethanol + 1min LABORETTE 17; SVM measurement

Calculation Fraunhofer broad **TradeOff** 1.000000e+002
Scans Fine 100 **Scans Coarse** 0 **Channels** 51 **Beam Obscuration** 15.00 %
Meas. Range 0.08 [µm] - 42.30 [µm] **Pump** 0 **Ultrasonics** 0

%	< µm
5	0.09
10	0.09
15	0.09
20	0.09
25	0.09
30	0.09
40	0.10
50	0.10
60	0.11
70	0.12
80	0.13
90	0.16
95	0.18
98	0.21
99	0.23

%	< µm
28.9	0.08
45.3	0.10
87.8	0.15
97.3	0.20
99.5	0.25
99.9	0.30
100.0	0.35
100.0	0.40
100.0	0.45
100.0	0.50
100.0	0.60
100.0	0.70
100.0	0.80
100.0	0.90
100.0	1.00

